

passion to lead
EPSILON

EP Corsa HF Mobile Solution



KW : 2.5 / 3.5 / 6 / 15 / 32

Certified & approved By:



Technical Specifications

SPECIFICATIONS	HIGH FREQUENCY MOBILE X-RAY GENERATOR					
MODEL	EP CORSA 2.5 M	EP CORSA 3.5 M	EP CORSA 6 SM	EP CORSA 6 M	EP CORSA 15 M	EP CORSA 30 M
Max Output (KW)	2.5	3.5	6	6	15	32
kV Range	40-110	40-110	40-110	40-125	40-125	40-125
mA Max.	60	80	150	150	300	400
X-Ray Tube Head	Monoblock Version					
X-Ray Tube	Stationary Anode			Rotating Anode		
Digital Display	KV, mA & mAs				KV & mAs	
Control Panel	Feather Touch (Micro Controller Based)					
APR Facility	Available					
Collimator	Manually / Motorized Light Beam Diaphragm (LBD)					
Power Required	1Ø 190-240 VAC 50 Hz with 15 Amps Power Socket (Required Proper Line Earthing)					
Tube Stands	Spring Balance Mobile / Counter Balance Mobile					
Battery/Capacitor Bank	Optional					

Due to continuous R & D for product development and improvement design, specifications, appearance & Data indicated might be vary from the product illustrated in the catalogue.



Marketed by: